Catchem'

POLYP TRAP

FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

EN - INSTRUCTIONS FOR USE
BG - ИНСТРУКЦИИ ЗА УПОТРЕБА
CS - NÁVOD K POUŽITÍ
DA - BRUGSANVISNING
DE - GEBRAUCHSANWEISUNG
EL - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ES - INSTRUCCIONES DE USO
ET - KASUTUSJUHISED
FI - KÄYTTÖOHJEET
FR - NOTICE
HR - UPUTE ZA UPORABU
HU - HASZNÁLATI UTASÍTÁS
IT - ISTRUZIONI PER L'USO
KK - ПАЙДАЛАНУ НҮСҚАУЛЫҒЫ
KO - 사용 지침
LT - NAUDOJIMO INSTRUKCIJOS
LV - LIETOŠANAS INSTRUKCIJA
NL - GEBRUIKSAANWIJZING
NO - BRUKSANVISNING
PL - INSTRUKCJA OBSŁUGI
PT - INSTRUÇÕES DE UTILIZAÇÃO
RO - INSTRUCȚIUNI DE UTILIZARE
RU - ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ
SK - NÁVOD NA POUŽITIE
SL - NAVODILA ZA UPORABO
SR - UPUTSTVO ZA UPOTREBU
SV - BRUKSANVISNING
UK - ІНСТРУКЦІЯ З ВИКОРИСТАННЯ


Endo-Therapeutics, Inc.
15251 Roosevelt Blvd. Suite 204
Clearwater, FL 33760 USA

EndoTherapeutics.com
Toll Free: **+1.844.400.ENDO** (3636)

Advena Ltd.
Tower Business Centre,
2nd Flr., Tower Street
Swatar, BKR 4013 Malta

Think Within.

EN – English.....	3
BG – Български (Bulgarian).....	3
CS – Čeština (Czech).....	4
DA – Dansk (Danish).....	4
DE – Deutsch (German).....	5
EL – Ελληνικά (Greek).....	5
ES – Español (Spanish).....	6
ET – Eesti (Estonian).....	7
FI – Suomi (Finnish).....	7
FR – Français (French).....	8
HR – Hrvatski (Croatian).....	8
HU – Magyar (Hungarian).....	9
IT – Italiano (Italian).....	9
KK – Қазақша (Kazakh).....	10
KO – 한국어 (Korean).....	11
LT – Lietuvių (Lithuanian).....	11
LV – Latviešu (Latvian).....	12
NL – Nederlands (Dutch).....	12
NO – Norsk (Norwegian).....	13
PL – Polski (Polish).....	13
PT – Português (Portuguese).....	14
RO – Română (Romanian).....	15
RU – русский (Russian).....	15
SK – Slovák (Slovak).....	16
SL – Slovenščina (Slovenian).....	16
SR – Srpski (Serbian).....	17
SV – Svenska (Swedish).....	17
UK – українська (Ukrainian).....	18
Symbols.....	19

EN – English

^{REF} **ET2211 Catchem'** Polyp Trap, 4 Chamber, Removable

^{REF} **ET2213 Catchem'** Polyp Trap, 4 Chamber, Fixed

Intended Use

The single use polyp trap is designed for in-line suction retrieval of polyps during endoscopic procedures.

Contraindications

None known within the scope of the intended use.

Warnings

Do not use if package is opened or damaged.

Visually inspect before use. If any abnormality is detected that would prohibit proper working condition, do not use.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure and/or transmission of disease.

Instructions

1. Connect suction tubing to center port on trap and outside tubing on lid to endoscope.
2. Rotate arrow on lid to position tubing over chamber 1.
3. Retrieve and collect specimen.
4. Rotate lid clockwise to position tubing over chamber 2 after first specimen is collected.
5. Retrieve and collect specimen.
6. Repeat step 4 and 5 by positioning tubing over chamber 3 & 4.
7. To detach, disconnect suction tubing from suction system and trap tubing from endoscope.
8. After all specimens are transferred to a transport vial, discard remaining portion of trap immediately.

Notes

Use of this device restricted to a trained healthcare professional.

In the event of a malfunction, discontinue use and save for return and contact your local representative.

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

Disposal

After use, this device may be a potential biohazard. Dispose of it in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Storage

Store in a dry location, away from temperature extremes.

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EndoT®, Think Within.® and Catchem® are registered trademarks of Endo-Therapeutics, Inc.
©2015-2021 Endo-Therapeutics, Inc. All rights reserved. 160403, Rev. F, 05/27/2021