

EndoTBasket

RETRIEVAL DEVICE

EndoTGrasper

RETRIEVAL DEVICE




STERILE EO



Rx only

FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

EN - INSTRUCTIONS FOR USE



Endo-Therapeutics, Inc.
15251 Roosevelt Blvd. Suite 204
Clearwater, FL 33760 USA

EndoTherapeutics.com
Toll Free: +1.844.400.ENDO (3636)

Think Within.



EN – English

<small>REF</small>	ET2030 EndoTBasket Elliptical Retrieval Basket, Medium
<small>REF</small>	ET2031 EndoTBasket Elliptical Retrieval Basket, Large
<small>REF</small>	ET2032 EndoTBasket Elliptical Retrieval Basket, Medium
<small>REF</small>	ET2033 EndoTBasket Elliptical Retrieval Basket, Large
<small>REF</small>	ET2034 EndoTBasket Elliptical Retrieval Basket, Medium
<small>REF</small>	ET2035 EndoTBasket Elliptical Retrieval Basket, Large
<small>REF</small>	ET2040 EndoTBasket Hexagonal Retrieval Basket, Medium
<small>REF</small>	ET2041 EndoTBasket Hexagonal Retrieval Basket, Large
<small>REF</small>	ET2042 EndoTBasket Hexagonal Retrieval Basket, Medium
<small>REF</small>	ET2043 EndoTBasket Hexagonal Retrieval Basket, Large
<small>REF</small>	ET2044 EndoTBasket Hexagonal Retrieval Basket, Medium
<small>REF</small>	ET2045 EndoTBasket Hexagonal Retrieval Basket, Large
<small>REF</small>	ET2020 EndoTGrasper Retrieval Device
<small>REF</small>	ET2021 EndoTGrasper Retrieval Device
<small>REF</small>	ET2022 EndoTGrasper Retrieval Device

Intended Use

Retrieval baskets and graspers are intended for use in removing foreign debris.

Contraindications

None known within the scope of the intended use.

Warnings and Precautions

This device is sharp. Use same care as for any other surgical device to avoid punctures or cuts to hands, gloves, etc.

Do not use if the sterile packaging is compromised, if any abnormalities are found or if the instrument is not operating properly.

Verify device is fully retracted into outer sheath before pulling device through endoscope channel.

Instructions

1. Peel open the pouch and remove the device. Visually inspect the device for kinks, loose or broken parts.
2. Form the device into three (3) 8-inch diameter loops for 220 cm length devices or two (2) 8 inch loops for 150cm length devices. Deploy and retract device several times to verify proper actuation.
3. Retract the device into the outer sheath before inserting the device into the working channel of the endoscope. If the endoscope is equipped with an elevator, raise the elevator fully before loading.
4. If resistance is met when the bending tip of the endoscope is angulated, release the bending angle slightly until the instrument passes smoothly. Do not force the device into an angulated endoscope as this could damage the endoscope.
5. After foreign debris is collected, retract device into outer sheath and slowly remove device from endoscope. If the scope has an elevator, it must be lowered before withdrawing the device.

Notes

Use of this device restricted to a trained healthcare professional.

In the event of a malfunction, discontinue use and save for return and contact your local representative.

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

Disposal

Dispose of the device after use in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Storage



Store sterile instruments at room temperature.

Do not expose to disinfectants, solvents, radiation, UV light or any other substance or environment that could compromise the sterility of the unit.

Do not place any objects on the device or packaging.

Symbols

EN - Symbols

	EN - Quantity of devices in pouch
	EN - Quantity of devices in box

Think Within.



Endo-Therapeutics, Inc.
15251 Roosevelt Blvd, Suite 204
Clearwater, FL 33760
USA

EndoTherapeutics.com
Toll Free: **+1.844.400.ENDO** (3636)

EndoT® and Think Within.® are registered trademarks of Endo-Therapeutics, Inc.
©2018-2024 Endo-Therapeutics, Inc. All rights reserved. P/N 151205, Rev F. 05/07/2024