eSuction

ENDOSCOPIC RETRIEVAL DEVICE

FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN

EN - INSTRUCTIONS FOR USE
BG - ИНСТРУКЦИИ ЗА УПОТРЕБА
CS - NÁVOD K POUŽÍTÍ
DA - BRUGSANVISNING
DE - GEBRAUCHSANWEISUNG
EL - ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
ES - INSTRUCCIONES DE USO
ET - KASUTUSJUHISED
FI - KÄYTTÖOHJEET
FR - NOTICE
HR - UPUTE ZA UPORABU
HU - HASZNÁLATI UTASÍTÁS
IT - ISTRUZIONI PER L'USO
KK - ПАЙДАЛАНУ НҮСҚАУЛЫҒЫ
KO - 사용 지침
LT - NAUDOJIMO INSTRUKCIJOS
LV - LIETOŠANAS INSTRUKCIJA
NL - GEBRUIKSAANWIJZING
NO - BRUKSANVISNING
PL - INSTRUKCJA OBSŁUGI
PT - INSTRUÇÕES DE UTILIZAÇÃO
RO - INSTRUCȚIUNI DE UTILIZARE
RU - ИНСТРУКЦИИ ПО ПРИМЕНЕНИЮ
SK - NÁVOD NA POUŽITIE
SL - NAVODILA ZA UPORABO
SR - UPUTSTVO ZA UPOTREBU
SV - BRUKSANVISNING
UK - ІНСТРУКЦІЯ З ВИКОРИСТАННЯ

Endo-Therapeutics, Inc.
15251 Roosevelt Blvd. Suite 204
Clearwater, FL 33760 USA

EndoTherapeutics.com
Toll Free: **+1.844.400.ENDO** (3636)

Advena Ltd.
Tower Business Centre,
2nd Flr., Tower Street
Swatar, BKR 4013 Malta

Think Within.

EN – English.....	3
BG – Български (Bulgarian).....	4
CS – Čeština (Czech).....	5
DA – Dansk (Danish).....	6
DE – Deutsch (German).....	7
EL – Ελληνικά (Greek).....	8
ES – Español (Spanish).....	9
ET – Eesti (Estonian).....	10
FI – Suomi (Finnish).....	11
FR – Français (French).....	12
HR – Hrvatski (Croatian).....	13
HU – Magyar (Hungarian).....	14
IT – Italiano (Italian).....	15
KK – Қазақша (Kazakh).....	16
KO – 한국어 (Korean).....	17
LT – Lietuvių (Lithuanian).....	18
LV – Latviešu (Latvian).....	19
NL – Nederlands (Dutch).....	20
NO – Norsk (Norwegian).....	21
PL – Polski (Polish).....	22
PT – Português (Portuguese).....	23
RO – Română (Romanian).....	24
RU – русский (Russian).....	25
SK – Slovák (Slovak).....	26
SL – Slovenščina (Slovenian).....	27
SR – Srpski (Serbian).....	28
SV – Svenska (Swedish).....	29
UK – українська (Ukrainian).....	30
Symbols.....	31

REF ET2000 eSuction Endoscopic Retrieval Device**REF ET2001 eSuction Endoscopic Retrieval Device****Indications**

The eSuction is intended to be used in the endoscopic retrieval of food bolus, foreign bodies and excised tissue such as polyps.

Contraindications

None known within the scope of the intended use.

Warnings and Precautions

Do not use if any abnormalities are found or if the instrument is not operating properly.

Verify device is fully retracted into outer sheath before pulling device through scope channel.

During retrieval, allow for proper spacing between the end of the eSuction cap and the object being retrieved before deploying the snare. Deploying the snare too closely to an object may result in the snare detaching from the anchor side of the eSuction cap.

This device is designed for single use only. Attempts to reprocess, resterilize, and/or reuse may lead to device failure, infection, and/or transmission of disease.

Do not use if the sterile packaging is compromised, any abnormalities are found or if the instrument is not operating properly.

Instructions

1. Verify that the outer diameter of the endoscope is compatible with the inner diameter of the eSuction cap. (Refer to device label.)
2. Peel open the pouch and remove the device. Visually inspect the device for kinks, loose or broken parts.
3. Deploy and retract the snare to ensure the device is functioning properly.
4. Verify that the endoscope has a 2.8 mm accessory channel or larger.
5. Dry the tip of the endoscope with a 4x4 gauze or dry wash cloth. It is important that the eSuction cap and the tip of the endoscope are dry to ensure that the device stays attached. Do not use any lubricant; doing so could result in the device falling off the endoscope inside the patient.
6. Load the retrieval snare into the accessory channel of the endoscope.
7. Position the distal tip of the snare into the catheter channel on the eSuction cap so that the snare distal tip is flush with the cap's distal end.
8. Place the eSuction cap onto the endoscope and align the catheter channel with the endoscope's accessory channel.
9. Gently pull the eSuction cap to ensure there is a proper fit on the endoscope.
10. Deploy the snare, fold over and attach the tip onto the anchor side of the eSuction cap.
11. Retract and deploy the snare to ensure attachment to anchor side of the eSuction cap.
12. The snare should be retracted before inserting into patient
13. Deploy snare loop around food bolus, foreign body or excised tissue by moving finger loops forward.
14. Retract snare loop with food bolus, foreign body or excised tissue by moving finger loops backward.
15. Endoscopic suction may be used as necessary to aid in the removal of food bolus impactions, foreign bodies and excised tissue.
16. After the procedure is completed, remove the endoscope from the patient.
17. Deploy and detach the snare loop from the anchor side of the eSuction cap and retract the snare.
18. Remove the eSuction cap from the endoscope and discard.

Notes

Use of this device restricted to a trained healthcare professional.

In the event of a malfunction, discontinue use and save for return and contact your local representative.

Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where the incident occurred.

Disposal

After use, this device may be a potential biohazard. Dispose of it in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Storage

Store sterile instruments at room temperature. Do not expose to solvents, radiation, UV light or any other any other substance or environment that could compromise the sterility of the unit.

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EndoT® and Think Within.® are registered trademarks of Endo-Therapeutics, Inc.
©2026 Endo-Therapeutics, Inc. All rights reserved. U.S. Patent 9,549,761 – 151009, Rev F