


# eSuction Endoscopic Retrieval Device



eSuction combines snare technology with suction to aid in the retrieval of food bolus impactions, foreign bodies and excised tissue such as polyps. The eSuction cap is designed to fit securely onto the endoscope.

- Food Bolus Impactions
- Foreign Bodies
- Excised Tissue

## Features and Benefits

---

- Snare Technology
  - Suction Accessory Cap
  - Sized for use in Upper and Lower GI procedures
- 

Think Within.




Endo-Therapeutics, Inc.  
15201 Roosevelt Blvd, Suite 104  
Clearwater, FL 33760

[EndoTherapeutics.com](http://EndoTherapeutics.com)

Toll Free: **844.400.ENDO** (3636)

# eSuction

## Endoscopic Suction Assisted Retrieval of Foreign Bodies



### ESOPHAGUS

#### Food Bolus Impactions


- Compressed meat
- Bone



### STOMACH

#### Foreign Bodies

- Ingested foreign objects
- Sharps
- Long object retrieval



### COLON

#### Excised Tissue

- Polyps
- Large lesions

**ET2000** eSuction, Esophageal - Fits Scope (O.D.) 8.6 mm – 10.0 mm - (3/box)

*For use with endoscope having a 2.8 mm accessory channel or larger.*

**ET2001** eSuction, Colonic - Fits Scope (O.D.) 10.8 mm – 13.2 mm - (3/box)

*For use with endoscope having a 2.8 mm accessory channel or larger.*



EndoTherapeutics.com • Toll Free: **844.400.ENDO** (3636)  
Endo-Therapeutics, Inc. • 15201 Roosevelt Blvd, Suite 104 • Clearwater, FL 33760

EndoT<sup>®</sup>, Think Within.<sup>®</sup> and eSuction<sup>®</sup> are registered trademarks of Endo-Therapeutics, Inc. U.S. Patent 9,549,761 ©2015-2018 Endo-Therapeutics, Inc. All rights reserved. MKT-008 Rev C 0118