

SwirlNet

Endoscopic Retrieval Device

SwirlNet is an endoscopic retrieval net that offers a safe and effective way of retrieving foreign bodies.

- Food Bolus Impactions
- Foreign Bodies
- Excised Polyp Removal

Features and Benefits

- Multiple sizes to accommodate various retrieval needs
 - Allows for easier removal of challenging foreign objects
 - Capable of handling multiple foreign body removals
 - Rotation for maximum control
-

Think Within.



Endo-Therapeutics, Inc.
15251 Roosevelt Blvd, Suite 204
Clearwater, FL 33760

EndoTherapeutics.com

Toll Free: **844.400.ENDO** (3636)

SwirlNet


Designed with soft flexible mesh for safe removal of foreign bodies, piecemeal polypectomy or endoscopic mucosal resection.



ESOPHAGUS

Food Bolus Impactions

- Compressed meat
- Bone



STOMACH

Foreign Bodies

- Ingested foreign objects
- Marbles, coins and dice
- Buttons and batteries



COLON

Excised Tissue

- Polyps
- Large lesions

Number	Description	Loop Size	Min. Channel	Tip Diameter	Length	Units
ET2010	SwirlNet Mini, Esophageal	2 x 4.5 Oval	2.0 mm	1.8 mm	160 cm	5
ET2011	SwirlNet Standard, Esophageal	3 x 6 Oval	2.8 mm	2.5 mm	160 cm	5
ET2012	SwirlNet Standard, Colonic	3 x 6 Oval	2.8 mm	2.5 mm	230 cm	10
ET2013	SwirlNet Standard, Enteroscope	3 x 6 Oval	2.8 mm	2.5 mm	350 cm	5
ET2014	SwirlNet Maxi, Esophageal	4 x 8 Oval	3.2 mm	3.0 mm	160 cm	5
ET2015	SwirlNet Maxi, Colonic	4 x 8 Oval	3.2 mm	3.0 mm	230 cm	5

All products are Non-Sterile



EndoTherapeutics.com • Toll Free: **844.400.ENDO** (3636)
Endo-Therapeutics, Inc. • 15251 Roosevelt Blvd, Suite 204 • Clearwater, FL 33760

EndoT[®], Think Within[®] and SwirlNet[®] are registered trademarks of Endo-Therapeutics, Inc. ©2017-2019 Endo-Therapeutics, Inc. All rights reserved. MKT-011 Rev C