

EndoT Needle

Performance matters.

Don't get stuck during your next injection therapy case



**Single Hand
Thumb Actuated
Needle Deployment**



**Triple Grind
Lancet Point**

**Stainless Steel
Sheath**

Think Within.



Endo-Therapeutics, Inc.
15251 Roosevelt Blvd, Suite 204
Clearwater, FL 33760


EndoTherapeutics.com

Toll Free: **844.400.ENDO** (3636)

EndoTNeedle

Use for tattooing, esophageal and colonic varices, ulcers and saline assisted polypectomy.

Designed for safety and ease of use, the spring return ergonomic handle offers a reliable and stable needle projection along with security from sharp injuries and scope damage.



Features and Benefits

Features a smooth round tip which aids in the closure and blockage needed to stop challenging bleed cases.

- Sterile
- Single-use, disposable
- Made in the USA

Number	Description	Needle Length	Min. Channel	Tip Diameter	Length	Units
ET1300	EndoTNeedle Endoscopic Injection Needle, 25 gauge	4 mm	2.8 mm	2.5 mm	160 cm	10
ET1301	EndoTNeedle Endoscopic Injection Needle, 25 gauge	4 mm	2.8 mm	2.5 mm	230 cm	10
ET1302	EndoTNeedle Endoscopic Injection Needle, 25 gauge	5 mm	2.8 mm	2.5 mm	160 cm	10
ET1303	EndoTNeedle Endoscopic Injection Needle, 25 gauge	5 mm	2.8 mm	2.5 mm	230 cm	10



EndoTherapeutics.com • Toll Free: **844.400.ENDO** (3636)
Endo-Therapeutics, Inc. • 15251 Roosevelt Blvd, Suite 204 • Clearwater, FL 33760