

Catchem' POLYP TRAP



ET2211

Catchem' PolypTrap, 4 Chamber, Removable

ET2213

Catchem' PolypTrap, 4 Chamber, Fixed

Features and Benefits

Packaged individually

Single-use device minimizes the risk of cross contamination.

In-line installation

Conveniently connects between the endoscope and suction canister.

Clear plastic design

Allows for visual confirmation during polyp collection.

4-Chamber index

Accurately identify each specimen sample in fixed or removable chambers.

Think Within.



Endo-Therapeutics, Inc.
15251 Roosevelt Blvd, Suite 204
Clearwater, FL 33760

EndoTherapeutics.com

Toll Free: **844.400.ENDO** (3636)

Capture excised tissue with a quick and convenient polyp collection.

ET2211

Catchem' Polyp Trap, 4 Chamber, Removable

Quickly and easily collect tissue samples using the removable chamber design.



ET2213

Catchem' Polyp Trap, 4 Chamber, Fixed

Traditional tissue collection is available using the fixed chamber design.

Number	Description	Units
ET2211	Catchem' Polyp Trap, 4 Chamber, Removable	12
ET2213	Catchem' Polyp Trap, 4 Chamber, Fixed	12



EndoTherapeutics.com • Toll Free: **844.400.ENDO** (3636)
Endo-Therapeutics, Inc. • 15251 Roosevelt Blvd, Suite 204 • Clearwater, FL 33760

EndoT®, Think Within.® and Catchem® are registered trademarks of Endo-Therapeutics, Inc. ©2015-2025 Endo-Therapeutics, Inc. All rights reserved. MKT-007 Rev D