

SINGLE-USE FLEXIBLE

# RhinoLaryngoscope

Used to examine and treat the nasal cavity and upper respiratory tract



Compatible with the MS-8 Image Processor

Working Length	Deflection	Field of View
350 mm	210° up/down	120°



### Increased Safety

Single-use eliminates risk of cross-infection



### Intensive Bending

210° up and down bending angle for articulated view



### Always Ready to Go

For immediate use and no reprocessing wait time



### User-friendly

Self-locking bending function and smooth angle control



### Portable

Lightweight for ease of use and transportation



### Comfortable Experience

Soft insertion tube for enhanced patient comfort

Think Within.

SINGLE-USE FLEXIBLE

# RhinoLaryngoscope

## Ergonomic Handle

Lightweight handle (under 300g) designed for comfort and ease of use

## Working Channel Port

Suction and biopsies can be performed through the working channel

## Suction Button

## Multifunction Buttons

Dual buttons for capturing photo and video, freezing frame, and white balance adjustments

## Distal Tip

Contains camera and LED

## Suction Port

Allows for connection of suction tubing

## Bending Control

Smooth bending function enables control of the distal tip

## Soft Insertion Tube

Provides a balance of rigidity for structural support and flexibility for ease of navigation

Number	Description	Working Length	Deflection	Field of View	Units
<b>RH-M22</b>	Single-use Flexible RhinoLaryngoscope (RH)- 2.2 mm, No Channel	350 mm	210° up/down	120°	1
<b>RH-M32A</b>	Single-use Flexible RhinoLaryngoscope (RH)- 3.1 mm, No Channel	350 mm	210° up/down	120°	1
<b>RH-M32</b>	Single-use Flexible RhinoLaryngoscope (RH)- 3.1 mm, 1.2 mm	350 mm	210° up/down	120°	1
<b>RH-M40</b>	Single-use Flexible RhinoLaryngoscope (RH)- 4.0 mm, 2.0 mm	350 mm	210° up/down	120°	1
<b>RH-M40H</b>	Single-use Flexible RhinoLaryngoscope (RH)- 4.0 mm, 1.2 mm	350 mm	210° up/down	120°	1



## DISTRIBUTED BY:

EndoTherapeutics.com • Toll Free: **844.400.ENDO** (3636)

Endo-Therapeutics, Inc. • 15251 Roosevelt Blvd, Suite 204 • Clearwater, FL 33760

EndoT® and Think Within.® are trademarks of Endo-Therapeutics, Inc. ©2025 Endo-Therapeutics, Inc. All rights reserved. MKT-027 Rev A