

# Image Processor

Portable High-definition Touchscreen Display

## High-Quality Imaging


15.6" HD touchscreen with 1920x1080 px resolution

## External Device Connection

Used for data transfers, image output, and network connection

## Information Management

Create patient records and generate diagnostic reports



## Dual Scope Display

Provides real-time views for both flexible and blade-type scopes simultaneously



Think Within.



DISTRIBUTED BY:

Endo-Therapeutics, Inc.  
15251 Roosevelt Blvd, Suite 204  
Clearwater, FL 33760

EndoTherapeutics.com

Toll Free: **844.400.ENDO** (3636)

# Designed for portability and intuitive use



## Compact and Flexible Design

- Convenient to transport and set up
- Can be placed on a cart or table, or mounted on a stand
- Versatile for use in the OR, ICU, or OPD

Weight: 2.6kg, Dimensions: 386.4 (h) x 253.6 (w) x 38.5 (d) mm

## Patient Information Management

- Create patient records and generate diagnostic reports
- Use DICOM protocol to transfer recorded imaging to hospital PACS system



## Compatibility

- Compatible with single-use endoscopes distributed by EndoT
- Port for Gold Fingers and Type-C connectors
- USB3.0 for faster data transfer
- Wired or wireless network connection
- DVI, HDMI and SDI image output

Number	Description	Units
MS-8	Image Processor	1



### DISTRIBUTED BY:

EndoTherapeutics.com • Toll Free: **844.400.ENDO** (3636)  
Endo-Therapeutics, Inc. • 15251 Roosevelt Blvd, Suite 204 • Clearwater, FL 33760

EndoT® and Think Within® are trademarks of Endo-Therapeutics, Inc. ©2025 Endo-Therapeutics, Inc. All rights reserved. MKT-030 Rev A