

SINGLE-USE FLEXIBLE

# Bronchoscope

Used for Airway Examinations, Intubation, Medication Delivery, Secretion Management, Biopsies, Alveolar Lavage and Foreign Body Retrieval



Compatible with the MS-8 Image Processor

Working Length	Deflection	Field of View
600 mm	210° up/down	120°



### Increased Safety

Single-use eliminates risk of cross-infection



### Intensive Bending

210° up and down bending angle for articulated view



### Always Ready to Go

For immediate use and no reprocessing wait time



### User-friendly

Ergonomic handle design and smooth bending control



### Portable

Lightweight and compact; convenient for the ICU



### Cost-Effective

Improves cash flow; no maintenance or reprocessing

Think Within.

# More Sizes, More Possibilities

Number	Description	Units
BR-M22	Single-use Flexible Bronchoscope (BR)- 2.2 mm, No Channel	1
BR-M32	Single-use Flexible Bronchoscope (BR)- 3.1 mm, 1.2 mm	1
BR-M40H	Single-use Flexible Bronchoscope (BR)- 4.0 mm, 1.2 mm	1
BR-M40	Single-use Flexible Bronchoscope (BR)- 4.0 mm, 2.0 mm	1
BR-M52H	Single-use Flexible Bronchoscope (BR)- 5.0 mm, 2.2 mm	1
BR-M50	Single-use Flexible Bronchoscope (BR)- 5.0 mm, 2.8 mm	1
BR-M58	Single-use Flexible Bronchoscope (BR)- 5.8 mm, 3.4 mm	1
BR-M58H	Single-use Flexible Bronchoscope (BR)- 6.0 mm, 2.8 mm	1

## Slim Models


For intubating or examining petite patients, checking DL tube position, and reaching smaller bronchi

## Routine Models

Can pass through most adult ET tubes for routine intubation and bronchoscopy

## Therapeutic Models

Regular outer diameter scopes with large working channels for performing therapeutic procedures



## DISTRIBUTED BY:



EndoTherapeutics.com • Toll Free: **844.400.ENDO** (3636)  
Endo-Therapeutics, Inc. • 15251 Roosevelt Blvd, Suite 204 • Clearwater, FL 33760

EndoT® and Think Within® are trademarks of Endo-Therapeutics, Inc. ©2025 Endo-Therapeutics, Inc. All rights reserved. MKT-031 Rev A