



Ergonomic Handle

Multi-Band Ligator

For delivering clear, quick, and precise ligation



Visual Reminder for One Shot Left

Multiple Bands for Treating Multiple Varices in One Operation

Think Within.





Distributed by:

Endo-Therapeutics, Inc.
15251 Roosevelt Blvd, Suite 204
Clearwater, FL 33760


EndoTherapeutics.com

Toll Free: 844.400.ENDO (3636)



Multi-Band Ligator

Use in conjunction with an endoscope for ligation of esophageal varices or anorectal hemorrhoids.



Features and Benefits

The Multi-Band Ligator is designed to perform several ligations in one operation. The carrier is made of transparent material, which allows clear visibility for precise ligation. It features tactile and audible band release cues, and the black band is a visual reminder for when there is only one band left.

Number	Description	Trigger Cord O.D.	Endoscope O.D.	Units
BL18-4A	Multi-Band Ligator- 4 Bands	0.2 mm	8.6-13 mm	1
BL18-6A	Multi-Band Ligator- 6 Bands	0.2 mm	8.6-13 mm	1
BL18-7A	Multi-Band Ligator- 7 Bands	0.2 mm	8.6-13 mm	1
BL18-9A	Multi-Band Ligator- 9 Bands	0.2 mm	8.6-13 mm	1



Distributed by:

EndoTherapeutics.com • Toll Free: **844.400.ENDO** (3636)
Endo-Therapeutics, Inc. • 15251 Roosevelt Blvd, Suite 204 • Clearwater, FL 33760

EndoT® and Think Within.® are trademarks of Endo-Therapeutics, Inc. ©2025 Endo-Therapeutics, Inc. All rights reserved. MKT-033 Rev A