

# EndoTBasket

Single-use Endoscopic Retrieval Basket

**Atraumatic Tip**

**3-Wire Design**

**Ergonomic Handle**

**Rotatable**

Think Within.




Endo-Therapeutics, Inc.  
15251 Roosevelt Blvd, Suite 204  
Clearwater, FL 33760

[EndoTherapeutics.com](http://EndoTherapeutics.com)

Toll Free: **844.400.ENDO** (3636)

# EndoTBasket

Used for the removal of resected polyps,  
foreign bodies, and esophageal food impactions



## Features and Benefits

The **EndoTBasket** features a 3-wire basket design with an atraumatic tip and ergonomic handle. It offers rotational capability for enhanced manipulation and secure grasping of resected polyps, foreign bodies, and esophageal food impactions.

- Sterile
- Single-use, disposable
- Manufactured in the USA

Number	Description	Opening Width	Min. Channel	Tip Diameter	Working Length	Units
ET2046	EndoTBasket Elliptical Retrieval Basket, Small	1.3 cm	2.0 mm	1.8 mm	230 cm	1



EndoTherapeutics.com • Toll Free: **844.400.ENDO** (3636)  
Endo-Therapeutics, Inc. • 15251 Roosevelt Blvd, Suite 204 • Clearwater, FL 33760

EndoT® and Think Within.® are trademarks of Endo-Therapeutics, Inc. ©2025 Endo-Therapeutics, Inc. All rights reserved. MKT-034 Rev A