

EndoTExtractor

Endoscopic Retrieval Device
with Excellent Pull Strength
for Long Objects

**4 Strong Grasping Prongs
with Atraumatic Tips**

Conical Design

**Single-Hand Operation
with 360° Grasping Force**

Think Within.



Endo-Therapeutics, Inc.
15251 Roosevelt Blvd, Suite 204
Clearwater, FL 33760

EndoTherapeutics.com

Toll Free: **844.400.ENDO** (3636)

EndoTExtractor


Designed for long object removal, featuring a control wire and loop system that enables maximum grasping strength for secure removal of foreign bodies.



ESOPHAGUS

Food Bolus Impactions

- Compressed meat
- Bone



STOMACH

Foreign Bodies

- Long object retrieval
- Sharps
- Ingested foreign bodies



COLON

Excised Tissue

- Polyps
- Large lesions

Number	Description	Loop Width	Min. Channel	Tip Diameter	Length	Units
ET2070	EndoTExtractor	25 mm	2.8 mm	2.6 mm	230 cm	3