

SINGLE-USE FLEXIBLE

RhinoLaryngoscope

Used to examine and treat the nasal cavity and upper respiratory tract



Compatible with the MS-8 Image Processor

Working Length	Deflection	Field of View
350 mm	210° up/down	120°



Increased Safety

Single-use eliminates risk of cross-infection



Intensive Bending

210° up and down bending angle for articulated view



Always Ready to Go

For immediate use and no reprocessing wait time



User-friendly

Self-locking bending function and smooth angle control



Portable

Lightweight for ease of use and transportation



Comfortable Experience

Soft insertion tube for enhanced patient comfort

Think Within.

SINGLE-USE FLEXIBLE

RhinoLaryngoscope

Ergonomic Handle

Lightweight handle (under 300g) with rotation handle for precise maneuverability.

Working Channel Port

Suction and biopsies can be performed through the working channel

Suction Button

Multifunction Buttons

Dual buttons for capturing photo and video, freezing frame, and white balance adjustments

Distal Tip

Contains camera and LED

Suction Port

Allows for connection of suction tubing

ET Tube Securement Connector

Bending Control

Smooth bending function enables control of the distal tip

Soft Insertion Tube

Provides a balance of rigidity for structural support and flexibility for ease of navigation

Number	Description	Units
RH-M22	Single-use Flexible RhinoLaryngoscope - 2.2 mm, No Channel	5
RH-M32A	Single-use Flexible RhinoLaryngoscope (HD) - 3.1 mm, No Channel	5
RH-M40	Single-use Flexible RhinoLaryngoscope - 4.0 mm, 2.0 mm	5
RH-M32	Single-use Flexible RhinoLaryngoscope - 3.1 mm, 1.2 mm	5
RH-M40H	Single-use Flexible RhinoLaryngoscope (HD) - 4.0 mm, 1.2 mm	5

The RH-M32A and RH-M40H have a higher definition view, allowing for clearer observation.

Diagnostic

Slim insertion tubes for examining nasal, pharyngeal, and laryngeal cavities

Therapeutic

Scopes with large working channels for suction, biopsies, and medication delivery in therapeutic procedures



DISTRIBUTED BY:

EndoTherapeutics.com • Toll Free: **844.400.ENDO** (3636)

Endo-Therapeutics, Inc. • 15251 Roosevelt Blvd, Suite 204 • Clearwater, FL 33760

EndoT® and Think Within.® are trademarks of Endo-Therapeutics, Inc. ©2025 Endo-Therapeutics, Inc. All rights reserved. MKT-027 Rev B