

Image Processor

Portable High-definition Touchscreen Display

High-Quality Imaging


15.6" HD touchscreen with
1920x1080 px resolution

External Device Connection

Used for data transfers, image
output, and network connection

Information Management

Create patient records and
generate diagnostic reports



Dual Scope Display

Provides real-time views for both
flexible and blade-type scopes
simultaneously



Think Within.



DISTRIBUTED BY:

Endo-Therapeutics, Inc.
15251 Roosevelt Blvd, Suite 204
Clearwater, FL 33760

EndoTherapeutics.com

Toll Free: **844.400.ENDO** (3636)

Designed for portability and intuitive use



Compact and Flexible Design

- Convenient to transport and set up
- Can be placed on a cart or table, or mounted on a stand
- Versatile for use in the OR, ICU, or OPD

Weight: 2.6kg, Dimensions: 386.4 (h) x 253.6 (w) x 38.5 (d) mm

Patient Information Management

- Create patient records and generate diagnostic reports
- Use DICOM protocol to transfer recorded imaging to hospital PACS system



Compatibility

- Compatible with single-use endoscopes distributed by EndoT
- Port for Gold Fingers and Type-C connectors
- USB3.0 for faster data transfer
- Wired or wireless network connection
- DVI, HDMI and SDI image output

Number	Description	Units
MS-8	Image Processor	1



DISTRIBUTED BY:

EndoTherapeutics.com • Toll Free: **844.400.ENDO** (3636)
Endo-Therapeutics, Inc. • 15251 Roosevelt Blvd, Suite 204 • Clearwater, FL 33760

EndoT® and Think Within® are trademarks of Endo-Therapeutics, Inc. ©2025 Endo-Therapeutics, Inc. All rights reserved. MKT-030 Rev B